




شرکت مهندسی پارس نهند
PARS NAHAND ENGINEERING Co.



Parsonic 2600s Ultrasonic Cleaner



Operating Voltage	220 V 50 Hz
Body, Tank, Basket & Drain Valve	Stainless Steel 304
Tank Capacity & Dimension	2.6 Lit. 240x135x100h mm
Body Dimension	330x165x210h mm
Ultrasonic Power & Frequency	70 Watt Max , 28 Khz ±5%
Heater	100 Watt 20-70 °C , Step 10 °C
Timer	2-22 min , Step 4 min
Accessory	Cable, Lid, Basket, Drain Hose
Guarantee	One Year

ISO 9001:2015 & 13485:2016 Certificate

Iran Ministry of Health Certificate for Medical Application

تهران ، کیلومتر ۲۰ جاده آبعلی ، پارک فن آوری پردیس ، خیابان نوآوری پنجم
پلاک ۵۱ کدپستی: ۱۶۵۴۱۲۰۸۷۷ تلفن/ فاکس: ۵ - ۰۲۱ - ۷۶۲۵۰۱۸۱
No. 51 , 5th Noavari Str. , Pardis TechPark , Tehran , Iran
Zip Code: 1654120877 TEL / FAX: +98 21 76250181-5
www.pnec.co.ir info@pnec.co.ir

