




شرکت مهندسی پارس نهند
PARS NAHAND ENGINEERING Co.



Parsonic 60sb Ultrasonic Cleaner



Operating Voltage	220 V 50 Hz
Body, Basket / Tank, Drain Valve	Stainless Steel 304 / 316L
Tank Capacity & Dimension	63 Lit. 700x300x300h mm
Body Dimension	895x385x540h mm
Ultrasonic Power & Frequency	600 Watt Max , 28 KHz ±5%
Heater	600 Watt 20-70 °C , Step 0.1 °C Digital
Timer	Digital Timer h/m/s
Accessory	Cable, Lid, Basket, Drain Hose
Guarantee	One Year

ISO 9001:2015 & 13485:2016 Certificate

Iran Ministry of Health Certificate for Medical Application

تهران ، کیلومتر ۲۰ جاده آبعلی ، پارک فن آوری پردیس ، خیابان نوآوری پنجم
پلاک ۵۱ کدپستی: ۱۶۵۴۱۲۰۸۷۷ تلفن/ فاکس: ۵ - ۷۶۲۵۰۱۸۱ - ۰۲۱
No. 51 , 5th Noavari Str. , Pardis TechPark , Tehran , Iran
Zip Code: 1654120877 TEL / FAX: +98 21 76250181-5
www.pnec.co.ir info@pnec.co.ir

